




extriCARE™ NPWT System - Comparison Chart

Differentiator	Smith & Nephew RENASYS™ GO	Smith & Nephew RENASYS™ EZ	KCI V.A.C. Freedom
			
Dimensions W x H x D	Oval Shape	14.5x9.5x7	6.5x3.125x7.5
Weight (lbs.)	2.4	7.4	3.2
Pressure Range	40~200mmHg	40~200mmHg	50~200mmHg
Canister	300cc & 750cc	250cc & 800cc	300cc
Vacuum Modes	Continuous or Intermittent	Continuous or Intermittent	Continuous or Intermittent
Battery Life	Lasts 20 hours Rechargeable	Lasts 40 hours Rechargeable	Lasts 12 hours Rechargeable
Major Alarms	<ul style="list-style-type: none"> • High/Low/Over Vacuum • Blockage/Canister Full • Low Battery • High Flow • Inactive 	<ul style="list-style-type: none"> • Over Vacuum • Blockage/Canister Full • Low Battery • High Flow/Leak 	<ul style="list-style-type: none"> • Blockage/Canister Full • Canister Installation • Low Battery • Leakage • Inactive
Unique Features	<ul style="list-style-type: none"> • Quiet and comfortable design • Easy transition from acute to home use • Intuitive operation 	<ul style="list-style-type: none"> • High flow rate • Long battery life • Multiple mounting options 	<ul style="list-style-type: none"> • On-screen user guide • TRAC technology simplifies dressing changes • Adjustable rate of exudates removal on dressing

extriCARE™ NPWT System - Comparison Chart

Differentiator	Devon Medical Products extriCARE™ 2400	KCI V.A.C. Via™	KCI ActiV.A.C.	Smith & Newpew PICO*
				
Dimensions W x H x D	1.46x5.67x3.35	6.8x3.9x2.4	7.6x6x2.5	3.5x3.5x1
Weight (lbs.)	.054	.07	2.4	0.27
Pressure Range	40~140mmHg	75 or 125mmHg	25~200mmHg	0~100mmHg
Canister	100cc, 250cc, 400cc	250cc	300cc	No Canister Used*
Vacuum Modes	Continuous or Intermittent	Continuous or Intermittent	Continuous or Intermittent	Continuous
Battery Life	Lasts 24 hours Rechargeable	Lasts 9 hours Rechargeable	Lasts 14 hours Rechargeable	AA 3V Battery
Major Alarms	<ul style="list-style-type: none"> • Air Leakage • Pressure Error • Canister Full • Tilt Sensor • Canister Installation • Low Battery 	<ul style="list-style-type: none"> • Blockage • Leakage • Low Battery • Therapy End 	<ul style="list-style-type: none"> • Blockage • Leakage • Low Battery • Canister Full 	<ul style="list-style-type: none"> • Leakage • Low Battery • Therapy End
Unique Features	<ul style="list-style-type: none"> • Easy & quick application • 5 Anatomically fitted dressings • Lightweight & portable • Long battery life 	<ul style="list-style-type: none"> • Uninterrupted therapy at discharge • Disposable after 7- day therapy 	<ul style="list-style-type: none"> • Color interface with touch-screen • Seal check for patient to self-trouble shoot • Patient therapy history access via USB port 	<ul style="list-style-type: none"> • Lightweight & portable • Disposable after 7-day therapy • No canister used*

* Not Medicare Reimbursable