


# An Advanced NPWT Dressing for Smaller Wounds



## ***The Same Kalypto Medical NPD 1000™ Technology***

Now available in a reduced size for smaller wounds with low drainage.

As with all NPD 1000 dressings, adjacent wounds can be treated simultaneously.

# Kalypto Medical NPD 1000™ NPWT System


## Portable. Effective. Easy.



### 25cc Dressing:

- Easy-to-Use, advanced dressing technology
  - "Peel and Stick" application
  - "Wicks and Locks" exudate
- Designed for smaller, low drainage wounds



### Dressing Size Comparison



25cc dressing      50cc dressing

### NPD 1000 Pump:

- Portable, lightweight and discreet
  - May increase patient compliance
  - Powered 3+ weeks by 3 AA batteries
  - No canister required

Product	Order #	Name	Description
	M1-05x08-S05	25cc Dressing; Case of 5	25cc Fluid Capacity. Dressing Size: 13 cm x16 cm Wound Contact Area: 5 cm x 8 cm. Pigtail Length: 6"
	P-1000i-K	NPD 1000 NPD 1000 System Device	NPD 1000 Pump & Control Module, Weight: 8 oz Dimensions: 3.2" W x 5.1" H x 1.3" D Device Tubing Length: 18" 3 AA Batteries

Contact your local Kalypto Medical Representative today for additional information on the NPD 1000 NPWT System.

### Kalypto Dressing Applications



25cc Dressing

50cc Oval

140cc Dressing

50cc Forefoot

70cc Sacral

